XM100/XM120 Platinum Electrodes

Storage

Replace the protection cap, and store the electrode in its box at room temperature.

Accessories

BAO3 Finishing Abrasive Bands	B90X550
ORENOVO•N	S16M001
Normal Cleaning Solution, 250 ml	A94L114
CL114 Electrode Cable, (1 m) with BNC plug	A94L116
CL116 Electrode Cable, (1 m) with type 7 plug	A94L116
CL111 Electrode cable, (1 m) with banana plug	A94L111

Specifications

Temperature range:	-10 to 100°C
Sensing element:	XM100: platinum wire ($\phi = 1 \text{ mm}$)
	XM120: platinum plate ($\phi = 5 \times 5 \text{ mm}$)

Introduction

The XM100 and XM120 Platinum Electrodes are ideal for a wide range of applications including redox measurements and redox titrations.

Preparation for Measurement

1. Remove the protection cap from the electrode.
2. Before starting a measurement, rinse the electrode with distilled water. Take care not to wet the screw cap.
3. Check that the electrode body and sensing element are not damaged.

Maintenance

Platinum rod or plate

The sensing part of the electrode can be polished with finishing abrasive bands, part no. B90X550. While polishing, take great care not to bend the platinum rod or plate.

In order to have reproducible and accurate results with the electrode, a few physical and electrochemical treatments are required. These treatments will depend on the last measurement performed:

Mild cleaning

Place the electrode overnight in RENOVO•N Normal Cleaning Solution.

Degreasing


- Remove mineral oils by rinsing the electrode with trichlorethylene and dabbing dry with a soft tissue.
- Remove organic oils by cathodic electrolysis in a hot alkaline solution, e.g. 15 g/l NaOH + 50 g/l Na_2CO_3 at 25°C.

Degassing

- Hydrogen produced during the electrolysis (see above) must be removed. Dip the electrode in water at 60°C and flush with pure nitrogen.

Chemical cleaning

- The electrode can be cleaned with a sulphochromic solution, or with a boiling solution of hydrochloric acid.

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D31M086 • Radiometer Analytical SAS • France • 9809A
RADIOMETER ANALYTICAL SAS